




# VALENCIA COLLEGE


## Chemistry

### Lab Technique 19: Using the LabQuest Data Collector with Conductivity Probe


Version 2

<p><b>Setup:</b></p> <p>Connect the conductivity probe into CH 1 port.</p> <p>Set the range switch on the amplifier box to the 0-20000 <math>\mu\text{S}</math> setting (the up position).</p> <p>Turn on the LabQuest Data Collector. Once it boots up it should automatically display CH 1: Conductivity (though the actual conductivity value will probably be different from 194 <math>\mu\text{S}/\text{cm}</math>).</p>	 <p>ON/OFF on the side</p> <p>Range switch up</p> <p>CH 1 port</p>
<p>While it boots up, suspend the conductivity probe from the stand using a utility clamp.</p> <p>Rinse the probe with deionized water (collect the rinse in the Waste beaker). Dry it well with a KimWipe, even inside of the elongated opening at the bottom (dispose of these in the regular garbage can).</p> <p>Once the conductivity probe is dry, leave it suspended on the clamp, in the air.</p>	
<p><b>Calibrate the Conductivity Probe<sup>1</sup>:</b></p>	
<p>Using the stylus connected to the LabQuest Data Collector, select Sensors from the tool bar on the top.</p>	 <p>Sensors</p>


Select Calibrate from the drop-down menu.



Click on CH 1: Conductivity 20000M|CS.




Click on Calibrate Now.



A new screen with a number pad on the lower portion will appear. Observe the reading displayed in the **Live Voltage** field. If the value is fluctuating, wait 10-15 seconds until it stabilizes.


With the dry conductivity probe suspended in the air, enter 0 in the box next to **Value 1** and Click the **Keep** box.




Lower the conductivity probe in the 1,000  $\mu\text{S}$  NaCl(aq) standard solution. (Leave this solution in the bottle provided and do NOT pour it into a beaker). Make sure that the probe's opening is submerged in the solution and swirl the bottle to remove any trapped air bubbles.

Once the voltage stabilizes, enter 1,000 in the box next to Known **Value 2** and click on the Keep box.

Click **OK**.



The original screen should now be displayed. The conductivity value should be fluctuating around 1,000  $\mu\text{S}/\text{cm}$ . (If not, something is wrong. Try re-calibrating it or get help.)



Remove the probe from the bottle and place the cap back on the bottle.

Rinse and dry the probe as before.



#### To Use:

Rinse the sensor with DI water and dry it well with a KimWipe.

Place sensor inside of the aqueous sample.

Make sure the opening on the bottom of the sensor is submerged in the solution. Stir it slightly to remove air bubbles.

Allow it to stabilize. If it doesn't quite stabilize, after 30 seconds, record the average reading.

When done with all readings, turn off, rinse and dry the sensor.

Reference:

1. Nord, R.S. *Conductivity and Chemical Reactions* [Online].

<https://www.emich.edu/chemistry/genchemlab/documents/5-conductivity.pdf> (accessed March 3, 2017).